



# Human Kalirin/RAC1 GEF/GTPase complex



## A Target Enabling Package (TEP)

<b>Gene ID / UniProt ID / EC</b>	KALRN: <a href="#">8997</a> / <a href="#">P97924</a> RAC1: <a href="#">5879</a> / <a href="#">P63000</a>
<b>Target Nominator</b>	SGC Internal Nomination
<b>SGC Authors</b>	Janine L. Gray, Carmen Jimenez Antunez, Tobias Krojer, Michael Fairhead, Nicola Burgess-Brown, Frank von Delft and Paul E. Brennan
<b>Collaborating Authors</b>	Roger Goody <sup>1</sup> , Romain Talon <sup>2</sup>
<b>Target PI</b>	Paul Brennan (SGC Oxford)
<b>Therapeutic Area(s)</b>	Neurological Disorders
<b>Disease Relevance</b>	KALRN has been linked to neurological disorders
<b>Date approved by TEP Evaluation Group</b>	June 10, 2019
<b>Document version</b>	Version 2
<b>Document version date</b>	October 2020
<b>Citation</b>	Janine L Gray, Carmen Jimenez Antunez, Tobias Krojer, Michael Fairhead, Romain Talon, Nicola Burgess-Brown, ... Paul E. Brennan. (2019). Human Kalirin/RAC1 GEF/GTPase Complex. [Data set]. Zenodo. <a href="https://doi.org/10.5281/zenodo.3266933">https://doi.org/10.5281/zenodo.3266933</a>
<b>Affiliations</b>	1. Max Planck Institute of Molecular Physiology, Dortmund 2. Diamond Light Source

We respectfully request that this document is cited using the DOI reference as given above if the content is used in your work.

### USEFUL LINKS



Open Targets



ChEMBL



[KALRN / RAC1](#)

[RAC1](#)

(Please note that the inclusion of links to external sites should not be taken as an endorsement of that site by the SGC in any way)

### SUMMARY OF PROJECT

We have cloned, expressed, purified and crystallised the first DH GEF domain of Kalirin (gene: [KALRN](#)) in complex with RAC1 (gene: [RAC1](#)) as part of a [programme](#) to explore a new way of targeting GTPases. Fragment screening and X-ray crystallography has identified binders within the orthosteric binding site. A nucleotide exchange assay has been developed and GEF activity established.

## SCIENTIFIC BACKGROUND

**RAC1** is a member of the RHO subfamily of the RAS small GTPases. These are guanine nucleotide dependent molecular switches fundamental to numerous cell signalling pathways and cellular processes. **KALRN** is a guanine nucleotide exchange factor (GEF), a type of regulatory protein which activates RAC1 by catalysing the exchange of GDP for GTP (1). Within full length Kalirin(2), there are two GEF domains, DH1 and DH2, which target two individual members of the RHO subfamily of small GTPases, proteins RAC1 and RHOA (2). Kalirin is found predominantly in the brain. The most abundant isoform is Kalirin-7, a truncation of full length Kalirin which contains a single RAC1 DH GEF domain. It was found that through activation of RAC1, Kalirin mediates dendritic spine formation and actin remodelling. Changes in KALRN expression have therefore been linked to a number of psychiatric and neurodegenerative disorders associated with dendritic spine pathology (2-4).




Fig. 1. Domain sequence of full length Kalirin.

Disease	Genetic association	Postmortem evidence	Molecular association	Animal models
Schizophrenia	GWAS (5,6)  Rare missense mutations (7)	↓ mRNA (8) ↓ kalirin-7 protein (9) ↑ kalirin-9 (10)	DISC1 (11,12) NRG1 /ErbB4 (13,14) 5-HT <sub>2A</sub> (15) NMDAR (12) PSD-95 (12) PAK	KALRN KO (16) Kal7 KO (17)
Alzheimer's disease	--	↓ mRNA (18) ↓ protein (18,19)	iNOS (20) PAK (21) EphB2 (22)	--
ADHD	GWAS (23)	--	Cadherins (24)	--
Addiction	--	--	--	Kal7 KO, addiction in mouse model (25,26)
Huntington's disease	--	--	HAP1 (27)	--
Parkinson's disease	--	--	Synphillin-1 (28)	--
Stress	--	--	--	Chronic restraint stress in mice (29)
Ischemic stroke	Case-control (30)	--	--	Mouse model of ischemia (31)

Table 1 Reproduced from Remmers et al., 2014 (2). The literature mentioned within can be found in the references.

Objectives of this TEP:

- Demonstrate that Kalirin/RAC1 can be expressed, purified and crystallised as recombinant proteins
- Identify nucleotide phosphate competitive chemical starting points for a Kalirin/RAC1 inhibitor
- Develop an assay to assess Kalirin/RAC1 activity


## RESULTS – THE TEP

### Purified proteins

Using *E. coli* as the expression host, we expressed and purified samples of the 1<sup>st</sup> DH GEF domain of KALRN (Kalirin (1)), the 1<sup>st</sup> DH and PH domains of KALRN (Kalirin (2)) (**Fig. 1**) and truncated (by 15 residues) full length RAC1 (RAC1).

### Structures


We have obtained 12 high resolution diffraction crystal structures of Kalirin(1)/RAC1 either as a holo-GDP bound structure (1.64 Å resolution, PDB ID [5O33](#), **Fig. 2**) or with fragments bound (1.9-2.3 Å resolution, PDB IDs [5QQD](#), [5QQE](#), [5QQF](#), [5QQG](#), [5QQH](#), [5QQI](#), [5QQJ](#), [5QQK](#), [5QQL](#), [5QQM](#), [5QQN](#)).











**Fig. 2** Structure of GDP-bound Kalirin(1)/RAC1 (Kalirin(1) green; RAC1 light blue; GDP, brown), PDB ID: [5O33](#)



















### Chemical Matter







A fragment soaking campaign using the XChem fragment library on GDP-bound crystals delivered 1 hit which displaced GDP. Docking and subsequent soaking of related compounds found a further 10 analogues.



**Fig. 3** Overlaid structures of fragment hits in the GDP binding site of Kalirin/RAC1 (Kalirin(1) green, RAC1 light blue, PDB IDs 5QQD, 5QQE, 5QQF, 5QQG, 5QQH, 5QQI, 5QQJ, 5QQK, 5QQL, 5QQM, 5QQN)

PDBID	Ligand	Binding Location	Binding Pocket	Resolution (Å)
5QQD	 <b>Z56880342</b> (Original hit)			1.91
5QQE	 <b>SG000012</b> (follow-up)			1.95
5QQF	 <b>SG000013</b> (follow-up)			2.26

5QQG	 <p>SG000055 (follow-up)</p>			2.23
5QQH	 <p>SG000059 (follow-up)</p>			2.09
5QQI	 <p>SG000070 (follow-up)</p>			2.08
5QQJ	 <p>SG000086 (follow-up)</p>			1.9
5QQK	 <p>SG000089 (follow-up)</p>			2.24
5QQL	 <p>SG000096 (follow-up)</p>			2.25


5QQM	 SG000098 (follow-up)			2.02
5QQN	 SG000112 (follow-up)			2.26

**Table 2** Chemical structures of the fragment hit and follow-up compounds

**IMPORTANT:** Please note that the existence of small molecules within this TEP indicates only that chemical matter might bind to the protein in potentially functionally relevant locations. The small molecule ligands are intended to be used as the basis for future chemistry optimisation to increase potency and selectivity and yield a chemical probe or lead series. As such, the molecules within this TEP should not be used as tools for functional studies of the protein, unless otherwise stated, as they are not sufficiently potent or well-characterised to be used in cellular studies.

### Nucleotide exchange activity assay

A nucleotide exchange assay was developed for Kalirin/RAC1 to measure the exchange rate of the GDP inside the pocket with GTP in solution (**Fig. 4**).



**Fig. 4** The nucleotide exchange assay. (1) Loading BODIPY-GDP into RAC1. (2) Fixing BODIPY-GDP into the pocket. (3) Nucleotide exchange reaction of bound BODIPY-GDP with GTP in solution.

This procedure measures the exchange reaction rate of the guanine nucleotide bound to the GTPase RAC1 with GTP in solution. It can be used to compare the exchange activity of RAC1 when it is either in complex with its GEF interacting partner Kalirin or alone. The assay makes use of a fluorescently labelled GDP analogue containing a boron-dipyrromethene moiety attached to position 2' or 3' (BODIPY-GDP, ex/em: 485/510 nm).

The fluorescence yield from BODIPY-GDP is enhanced when it is bound to RAC1 and decreases upon release. This can be used to measure the nucleotide exchange rate through changes in the fluorescence intensity.

The assay has three phases: **(1)** nucleotide replacement, **(2)** BODIPY-GDP fixing, and **(3)** nucleotide exchange (Fig. 4). During the nucleotide replacement, EDTA sequesters the  $Mg^{++}$  ion, essential for GDP binding, and GDP is displaced by labelled BODIPY-GDP. When equilibrium is reached, an excess of  $Mg^{++}$  is added to complex the remaining EDTA and fix the fluorescent nucleotide inside the pocket. Lastly, GTP is added to trigger the dissociation of BODIPY-GDP and RAC1, decreasing the fluorescence intensity signal (Fig. 5)




Fig. 5. The nucleotide exchange assay shows that SGC-expressed Kalirin is active as a GEF. Solid line represents initial rate of Kalirin (2)/RAC1.

#### **Future work**

- Develop chemical probes for Kalirin/RAC1 based on TEP fragment hits
- Express and screen other GEF/GTPase targets

#### **Collaborations**

- A collaboration with Paul Harrison (University of Oxford) is planned to test compounds binding to Kalirin/Rac1 in cellular models of schizophrenia.
- Nir London, received a grant from the Israeli Cancer Research Fund to develop RAC1 inhibitors in collaboration with the Oxford team (<https://www.icrffonline.org/scientists/nir-london-phd/>).

## **CONCLUSION**

We have generated protein, an assay, crystal structures and chemical matter that has been shown to bind to the guanine nucleotide binding pocket of the Kalirin/RAC1 complex.

## **TEP IMPACT**

Publications arising from this work:

- Gray, J.L., von Delft, F. and Brennan, P.E. [Targeting the Small GTPase Superfamily through their Regulatory Proteins](#). (2019) *Angew. Chem. Int. Ed.*

## **FUNDING INFORMATION**

The work at the SGC has been supported by Wellcome (grant 106169/ZZ14/Z), the Innovative Medicines Initiative Joint Undertaking (IMI JU) (grant 115766), Alzheimer's Research UK (grant ARUK-2018DDI-OX), National Institutes of Health (NIH) (grant 1RF1AG057443-01) and NDM Prize Graduate Studentship (Nuffield Department of Medicine, University of Oxford and Diamond Light Source Ltd, UK)



## ADDITIONAL INFORMATION

### Structure Files

PDB ID	Structure Details	Compound ID
5O33	Structure of GDP-bound Kalirin/RAC1	GT000001f
5QQD	Kalirin/RAC1 in complex with fragment	Z56880342
5QQE	Kalirin/RAC1 in complex with fragment	SG000012
5QQF	Kalirin/RAC1 in complex with fragment	SG000013
5QQG	Kalirin/RAC1 in complex with fragment	SG000055
5QQH	Kalirin/RAC1 in complex with fragment	SG000059
5QQI	Kalirin/RAC1 in complex with fragment	SG000070
5QQJ	Kalirin/RAC1 in complex with fragment	SG000086
5QQK	Kalirin/RAC1 in complex with fragment	SG000089
5QQM	Kalirin/RAC1 in complex with fragment	SG000096
5QQN	Kalirin/RAC1 in complex with fragment	SG000098
5QQO	Kalirin/RAC1 in complex with fragment	SG000112

### Materials and Methods

#### **Protein expression and purification**

##### ***Kalirin (1)***

**Vector:** pNIC28-10His

**Cell line:** BL21(DE3)-R3-pRARE2

**Tags and additions:** N-terminal, TEV protease cleavable decahistidine tag

Final protein sequence:

MHHHHHHHHHHSSGVDLGTENLYFQ\*SMRKKEFIMAELLQTEKAYVRDLHECLETYLWEMTSGVEEIPPGILNKEHIIFGNI  
QEIYDFHNNIFLKELEKYEQLPEDVGHCFVTWADKFQMYVTYCKNKPDSNQLILEHAGTFFDEIQRHGLANSISSYLIKPVQ  
RVTKYQLLLKELLTCCEEGKGGELKDGLEVMLSVPPKKANDAMHV

(underlined sequence contains vector encoded His-tag and TEV protease cleavage site\*)

##### **Protein Expression**

A 10 mL overnight culture grown in LB containing kanamycin (final concentration 50 µg/mL) and chloramphenicol (34 µg/mL) at 37 °C was used to inoculate 1 L of AIM-TB (ForMedium) with 0.01% Antifoam 204 and antibiotics in a 2.5 L baffled flask (Ultrayield, Thomson). Cells were grown for 4 hours at 37 °C 250 rpm before the temperature reduced to 25 °C 250 rpm and shaken overnight. Cells were harvested by centrifugation at 4,000 x g for 20 minutes at 4 °C and the pellet stored at -20 °C.

##### **Cell Lysis**

Lysis Buffer: 10 mM HEPES pH 7.5, 500 mM NaCl, 5% Glycerol, 0.5 mM TCEP, 0.5 mg/mL Lysozyme, 0.1 µg/mL Benzonase and protease inhibitors (Calbiochem EDTA-free Protease Inhibitor Cocktail Set III)

The cell pellet was resuspended in lysis buffer (20 mL per 5g of pellet), Triton X-100 to a final concentration of 2% added and the cells frozen overnight at -80 °C. The cells were thawed in a room temperature water bath, imidazole added to a final concentration of 20 mM and the cell debris was removed by centrifugation for 1 hour at 4000 x g.

##### **Column 1: His GraviTrap columns (3 x 1 mL volume of Ni-Sepharose 6 FF in a gravity-flow column)**

Wash Buffer: 10 mM HEPES pH 7.5, 500 mM NaCl, 5% Glycerol, 20 mM Imidazole, 0.5 mM TCEP

Elution Buffer: 10 mM HEPES pH 7.5, 500 mM NaCl, 5% Glycerol, 500 mM Imidazole, 0.5 mM TCEP

The clarified cell extract was added to 3 x 1 mL of His GraviTrap columns pre-equilibrated with Wash buffer. The columns were then washed with 2 x 10 mL Wash Buffer. The protein was eluted with 2.5 mL Elution Buffer after applying column 1 directly onto column 2.

Column 2: PD-10 desalting columns (3 x 8.3 mL volume of Sephadex G-25 resin in a gravity-flow column)

Wash Buffer: 10 mM HEPES pH 7.5, 500 mM NaCl, 5% Glycerol, 20 mM Imidazole, 0.5 mM TCEP

Following elution of the protein from the His GraviTrap column directly onto the PD-10 column, the HisGraviTrap column was removed and the protein eluted from the PD-10 column using 3.5 mL of Wash Buffer.

Tag cleavage

1 mg of TEV protease was added to every 10 mg of the eluted protein and the digestion was performed overnight at 4 °C.

Column 3: His GraviTrap columns (3 x 1 mL volume of Ni-Sepharose 6 FF in a gravity-flow column)

Wash Buffer: 10 mM HEPES pH 7.5, 500 mM NaCl, 5% Glycerol, 20 mM Imidazole, 0.5 mM TCEP

The TEV cleaved protein was applied to a HisGraviTrap column pre-equilibrated in Wash Buffer to remove TEV protease, His-tag and uncleaved protein. The columns were then washed with 2.5 mL of Wash Buffer for a final pool of 6 mL. Protein fractions were combined and concentrated to ~ 30 mg/mL using a 10 kDa MWCO concentrator.

Column 4: Yarra SEC 2000 300x21.2 mm column (300x21.2 mm column, 104 mL volume) (Gel filtration)

GF buffer: 10 mM HEPES pH 7.5, 500 mM NaCl, 5% glycerol, 0.5 mM TCEP

The protein was loaded onto the column pre-equilibrated in GF buffer, the peak corresponding to the target protein was taken, concentrated to 12 mg/ml and stored at -80 °C.

**Kalirin (2)**

**Vector:** pNIC28-Bsa4

**Cell line:** BL21(DE3)-R3-pRARE2

**Tags and additions:** N-terminal, TEV protease cleavable hexahistidine tag

Final protein sequence:

MHHHHHHSSGVDLGTENLYFQ\*SMDREVKLRDANHEVNEEKRKSARKKEFIMAELLQTEKAYVRDLHECLETYLWEMTSG  
VEEIPPILNKEHIIIFGNIQEYDFHNNIFLKELEKYEQLPEDVGHCFVTWADKFMVYTYCKNKPDSNQLILEHAGTFFDEIQQ  
RHGLANSISSYLIKPVQRVTKYQLLLKELLTCCEEGKELKDGLEVMLSVPKKANDAMHVSMLLEGFDENLDVQGELILQDAFQ  
VWDPKSLIRKGRERHLFLFEISLVFSKEIKDSSGHTKYVYKNKLLTSELGVTEHVEGDPCKFALWSGRTSPSSDNKTVLKASNIETK  
QEWIKNIREVIQERIIHLK GAL

(underlined sequence contains vector encoded His-tag and TEV protease cleavage site\*)

Protein Expression

A 10 mL overnight culture grown in LB containing kanamycin (final concentration 50 µg/mL) and chloramphenicol (34 µg/mL) at 37 °C was used to inoculate 1 L of AIM-TB (ForMedium) with 0.01% Antifoam 204 and antibiotics in a 2.5 L baffled flask (Ultrayield, Thomson). Cells were grown for 4 hours at 37 °C 250 rpm before the temperature reduced to 25 °C 250 rpm and shaken overnight. Cells were harvested by centrifugation at 4,000 x g for 20 minutes at 4 °C and the pellet stored at -20 °C.

Cell Lysis

Lysis Buffer: 10 mM HEPES pH 7.5, 500 mM NaCl, 5% Glycerol, 0.5 mM TCEP, 0.5 mg/ml Lysozyme, 0.1 µg/ml Benzonase and protease inhibitors (Calbiochem EDTA-free Protease Inhibitor Cocktail Set III)

Cell pellet was resuspended in lysis buffer (20 mL per 5g of pellet), Triton X-100 to a final concentration of 2% added and the cells frozen overnight at -80 °C. The cells were thawed in a room-temperature water bath, imidazole added to a final concentration of 20 mM and the cell debris was removed by centrifugation for 1 hour at 4000 x g.

Column 1: His GraviTrap columns (3 x 1 mL volume of Ni-Sepharose 6 FF in a gravity-flow column)

Wash Buffer: 10 mM HEPES pH 7.5, 500 mM NaCl, 5% Glycerol, 20 mM Imidazole, 0.5 mM TCEP

Elution Buffer: 10 mM HEPES pH 7.5, 500 mM NaCl, 5% Glycerol, 500 mM Imidazole, 0.5 mM TCEP

The clarified cell extract was added to 3 x 1 mL of His GraviTrap columns pre-equilibrated with Wash buffer. The columns were then washed with 2 x 10 mL Wash Buffer. The protein was eluted with 2.5 mL Elution Buffer after applying column 1 directly onto column 2.

Column 2: PD-10 desalting columns (3 x 8.3 mL volume of Sephadex G-25 resin in a gravity-flow column)

Wash Buffer: 10 mM HEPES pH 7.5, 500 mM NaCl, 5% Glycerol, 20 mM Imidazole, 0.5 mM TCEP

Following elution of the protein from the His GraviTrap column directly onto the PD-10 column, the HisGraviTrap column was removed and the protein eluted from the PD-10 column using 3.5 mL of Wash Buffer.

Tag cleavage

1 mg of TEV protease was added to every 10 mg of the eluted protein and the digestion was performed overnight at 4 °C.

Column 3: His GraviTrap columns (3 x 1 mL volume of Ni-Sepharose 6 FF in a gravity-flow column)

Wash Buffer: 10 mM HEPES pH 7.5, 500 mM NaCl, 5% Glycerol, 20 mM Imidazole, 0.5 mM TCEP

The TEV cleaved protein was applied to a HisGraviTrap column pre-equilibrated in Wash Buffer to remove TEV protease, His-tag and uncleaved protein. The columns were then washed with 2.5 mL of Wash Buffer for a final pool of 6 mL. Protein fractions were combined and concentrated to ~ 30 mg/mL using a 30 kDa MWCO concentrator.

Column 4: Yarra SEC 2000 300x21.2 mm column (300x21.2 mm column, 104 mL volume) (Gel filtration)

GF buffer: 10 mM HEPES pH 7.5, 500 mM NaCl, 5% glycerol, 0.5 mM TCEP

The protein was loaded onto the column pre-equilibrated in GF buffer, the peak corresponding to the target protein was taken, concentrated to 20 mg/ml and stored at -80 °C.

**RAC1**

**Vector:** pNIC-NStIIT

**Cell line:** BL21(DE3)-R3-pRARE2

**Tags and additions:** N-terminal, TEV protease cleavable Strep-tag II

Final protein sequence

HMSSGASWHPQFEKGGGSGGGSGGAWSHPQFEKGSVDLGTENLYFQ\*SMQAIKCVVVDGAVGKTCLLISYTTNAF  
PGEYIPTVFDNYSANVMVDGKPVNLGLWDTAGQEDYDRLRPLSYPTDVFSLVSPASFENVRKRWYPEVRHHCPNTPI  
ILVGTKLDRDDKDTIEKLEKLLTPITYPQGLAMAKEIGAVKYLECSALTQRGLKTVFDEAIRAVL

(underlined sequence contains vector encoded Strep-tag and TEV protease cleavage site\*)

Protein Expression

A 10 mL overnight culture grown in LB containing kanamycin (final concentration 50 µg/mL) and chloramphenicol (34 µg/mL) at 37 °C was used to inoculate 1 L of AIM-TB (ForMedium) with 0.01% Antifoam 204 and antibiotics in a 2.5 L baffled flask (Ultrayield, Thomson). Cells were grown for 4 hours at 37 °C 250 rpm

before the temperature reduced to 25 °C 250 rpm and shaken overnight. Cells were harvested by centrifugation at 4,000 x g for 20 minutes at 4 °C and the pellet stored at -20 °C.

#### Cell Lysis

Lysis Buffer: 10 mM HEPES pH 7.5, 500 mM NaCl, 5% Glycerol, 0.5 mM TCEP, 0.5 mg/ml Lysozyme, 0.1 µg/ml Benzonase and protease inhibitors (Calbiochem EDTA-free Protease Inhibitor Cocktail Set III)

Cell pellet was resuspended in lysis buffer (20 mL per 5g of pellet), Triton X-100 to a final concentration of 2% added and the cells frozen overnight at -80 °C. The cells were thawed in a room temperature water bath and the cell debris was removed by centrifugation for 1 hour at 4000 x g.

#### Column 1: Strep-Tactin columns (3 x 1 mL volume of Strep-Tactin XT resin in a gravity-flow column)

Wash Buffer: 10 mM HEPES pH 7.5, 500 mM NaCl, 5% Glycerol, 0.5 mM TCEP

Elution Buffer: 10 mM HEPES pH 7.5, 500 mM NaCl, 5% Glycerol, 50 mM D-Biotin, 0.5 mM TCEP

The clarified cell extract was added to 3 x 1 mL of His GraviTrap columns pre-equilibrated with Wash buffer. The columns were then washed with 2 x 10 mL Wash Buffer. 2 mL of Elution Buffer was applied to each column. After 10 minutes and 20 minutes, this was repeated to result in a 6 mL pool of eluted protein.

#### Tag cleavage

1 mg of TEV protease was added to every 10 mg of the eluted protein and the digestion was performed overnight at 4 °C.

#### Column 2: His GraviTrap columns (3 x 1 mL volume of Ni-Sepharose 6 FF in a gravity-flow column)

Wash Buffer: 10 mM HEPES pH 7.5, 500 mM NaCl, 5% Glycerol, 20 mM Imidazole, 0.5 mM TCEP

The TEV cleaved protein was applied to a HisGraviTrap column pre-equilibrated in Wash Buffer to remove TEV protease. The columns were then washed with 2.5 mL of Wash Buffer. Protein fractions were combined and concentrated to ~ 30 mg/mL using a 10 kDa MWCO concentrator.

#### Column 3: Yarra SEC 2000 300x21.2 mm column (300x21.2 mm column, 104 mL volume) (Gel filtration)

GF buffer: 10 mM HEPES pH 7.5, 500 mM NaCl, 5% glycerol, 0.5 mM TCEP

The protein was loaded onto the column pre-equilibrated in GF buffer, the peak corresponding to the target protein was taken and concentrated to 15 mg/ml and stored at -80 °C.

### ***Kalirin(1)/RAC1 GDP-free complex formation***

#### GDP removal

HPLC buffer: 100 mM Potassium phosphate pH 6.5, 10 mM Tetrabutyl phosphonium bromide, 5% methanol

Equipment: Agilent 1200 series HPLC. Reverse-phase C-18 5 µM (250 x 4.6 mm)

10 x Exchange Buffer: 2M ammonium phosphate, 10 µM zinc chloride

Kalirin (1) and RAC1 were mixed in a 1:1 molar ratio for an hour and concentrated to 60 mg/mL (138 µL) using a 30 kDa MWCO concentrator. 40 µL of 10 mM β,γ-Methyleneguanosine 5'-triphosphate sodium salt (GppCp, Sigma) and 4 Units of Alkaline phosphatase (Roche) were added. 20 µL of the 10x exchange buffer was added (final concentration 1 x) and the mixture incubated at 4 °C for 3 hours (Final concentration of Kalirin(1)/RAC1 is 1 mM). Degradation of GDP to GMP and guanosine was monitored by HPLC. Following complete degradation of GDP, 0.008 Units of snake venom phosphodiesterase I (Sigma) was added. The reaction was left overnight at 4 °C and the degradation of GppCp to GMP and guanosine monitored by HPLC. The protein was diluted to 30 mg/mL, flash frozen in liquid nitrogen and stored at -80 °C.

#### Column 1: Yarra SEC 2000 300x21.2 mm column (300x21.2 mm column, 104 mL volume)

GF buffer: 10 mM HEPES pH 7.5, 500 mM NaCl, 5% glycerol, 0.5 mM TCEP

The protein was thawed and loaded onto the column pre-equilibrated in GF buffer, the peak corresponding to the complex taken, concentrated to 10.4 mg/mL and stored at -80 °C.

### **Protein crystallisation**

Kalirin(1)/RAC1 was crystallised by mixing 75nL of 10.4 mg/mL protein in 10mM HEPES pH 7.5, 500mM NaCl, 5% Glycerol 0.5 mM TCEP with 75nL of reservoir solution containing 0.1M bis-tris pH 5.5, 24% PEG3350. Crystals appeared overnight from sitting drop plates at 20 °C. Kalirin(1)/RAC1 crystallised in space group P65 2 2 with unit cell dimensions of a=63 Å, b=63 Å, c=346 Å, corresponding to one Kalirin(1)/RAC1 molecule in the asymmetric unit. The crystals typically diffract between 1.7 and 2.3 Å.

### **Nucleotide exchange assay**

#### **Equipment**

PHERASTAR FSX Plate reader (BMG Labtech, Germany)

#### **Consumables**

Microplate, 384 well, PS, flat bottom, small volume, non-binding, black, 784900. (Grenier Bio-One, UK)

#### **Chemicals and reagents**

MgCl<sub>2</sub> (magnesium chloride) M8266. (Sigma-Aldrich, UK)

GTP (guanosine 5'-triphosphate disodium salt trihydrate) JBS-NU-1012-100. (Enzo Life Sciences, UK)

#### **Buffers**

Dilution buffer (DB): 20 mM TRIS HCl pH7.5, 50 mM NaCl, 1 mM MgCl<sub>2</sub>, 0.1% BSA, 1 mM DTT

Exchange buffer (EB): 20 mM TRIS HCl pH7.5, 50 mM NaCl, 2 mM EDTA, 0.1% BSA, 1 mM DTT, 0.75 μM BODIPY™ FL GDP

#### **Protein and peptide**

RAC1 c001 (RAC1), SGC: stored in aliquots at -80°C.

Kalirin c002 (Kalirin (2)), SGC: stored in aliquots at -80°C.

RAC1 was prepared in a stock solution at 50 μM in DB, Kalirin (2) was prepared in stock solution at 1.25 μM in DB. Three vials were prepared: 2 μL of RAC1, 2 μL of Kalirin (2) and 63.5 μL of EB were mixed into Vial A; 2 μL of RAC1, 2 μL of Kalirin (2) and 63.5 μL of EB were mixed into Vial B; 4 μL of DB and 63.5 μL of EB were mixed into Vial C. The vials were incubated at room temperature for 20 minutes, then 7.5 μL of 50 mM MgCl<sub>2</sub> in DB was added to each vial and mixed. To prepare the assay plate, 15 μL of Vial A was dispensed into row A columns 1-4, 15 μL of Vial B was dispensed into row B columns 1-4 and 15 μL of Vial C was dispensed into row C columns 1-4. The Pherastar FSX was used to dispense 5 μL of 1.8 mM GTP in DB into each well. The plate was then shaken for 10 seconds and FI signal was measured every 30 seconds for 75 minutes. FI 485/520, Gain 300. FI data was analysed in Prism 7 (Version 7.04). The enzyme rate was calculated as the slope obtained by fitting a straight line in a range from 0-5 minutes, units are ΔRFU/time.

### **References**

1. Wennerberg, K., Rossman, K. L., and Der, C. J. (2005) [The Ras superfamily at a glance](#). *J Cell Sci* **118**, 843-846
2. Remmers, C., Sweet, R. A., and Penzes, P. (2014) [Abnormal kalirin signaling in neuropsychiatric disorders](#). *Brain Res Bull* **103**, 29-38
3. Penzes, P., and Remmers, C. (2012) [Kalirin signaling: implications for synaptic pathology](#). *Mol Neurobiol* **45**, 109-118

4. Sommer, J. E., and Budreck, E. C. (2009) [Kalirin-7: linking spine plasticity and behavior](#). *J Neurosci* **29**, 5367-5369
5. Jean, P. S. (2008) Genes associated with schizophrenia identifies using a whole genome scan. in *US Patent Office*, Smithkline Beecham Corporation, Philadelphia, PA (US), United States
6. Ikeda, M., Aleksic, B., Kinoshita, Y., Okochi, T., Kawashima, K., Kushima, I., Ito, Y., Nakamura, Y., Kishi, T., Okumura, T., Fukuo, Y., Williams, H. J., Hamshere, M. L., Ivanov, D., Inada, T., Suzuki, M., Hashimoto, R., Ujike, H., Takeda, M., Craddock, N., Kaibuchi, K., Owen, M. J., Ozaki, N., O'Donovan, M. C., and Iwata, N. (2011) [Genome-wide association study of schizophrenia in a Japanese population](#). *Biol Psychiatry* **69**, 472-478
7. Kushima, I., Nakamura, Y., Aleksic, B., Ikeda, M., Ito, Y., Shiino, T., Okochi, T., Fukuo, Y., Ujike, H., Suzuki, M., Inada, T., Hashimoto, R., Takeda, M., Kaibuchi, K., Iwata, N., and Ozaki, N. (2012) [Resequencing and association analysis of the KALRN and EPHB1 genes and their contribution to schizophrenia susceptibility](#). *Schizophr Bull* **38**, 552-560
8. Hill, J. J., Hashimoto, T., and Lewis, D. A. (2006) [Molecular mechanisms contributing to dendritic spine alterations in the prefrontal cortex of subjects with schizophrenia](#). *Mol Psychiatry* **11**, 557-566
9. Rubio, M. D., Haroutunian, V., and Meador-Woodruff, J. H. (2012) [Abnormalities of the Duo/Ras-related C3 botulinum toxin substrate 1/p21-activated kinase 1 pathway drive myosin light chain phosphorylation in frontal cortex in schizophrenia](#). *Biol Psychiatry* **71**, 906-914
10. Deo, A. J., Cahill, M. E., Li, S., Goldszer, I., Henteleff, R., Vanleeuwen, J. E., Rafalovich, I., Gao, R., Stachowski, E. K., Sampson, A. R., Lewis, D. A., Penzes, P., and Sweet, R. A. (2012) [Increased expression of Kalirin-9 in the auditory cortex of schizophrenia subjects: its role in dendritic pathology](#). *Neurobiol Dis* **45**, 796-803
11. Millar, J. K., Christie, S., and Porteous, D. J. (2003) [Yeast two-hybrid screens implicate DISC1 in brain development and function](#). *Biochem Biophys Res Commun* **311**, 1019-1025
12. Hayashi-Takagi, A., Takaki, M., Graziane, N., Seshadri, S., Murdoch, H., Dunlop, A. J., Makino, Y., Seshadri, A. J., Ishizuka, K., Srivastava, D. P., Xie, Z., Baraban, J. M., Houslay, M. D., Tomoda, T., Brandon, N. J., Kamiya, A., Yan, Z., Penzes, P., and Sawa, A. (2010) [Disrupted-in-Schizophrenia 1 \(DISC1\) regulates spines of the glutamate synapse via Rac1](#). *Nat Neurosci* **13**, 327-332
13. Cahill, M. E., Jones, K. A., Rafalovich, I., Xie, Z., Barros, C. S., Muller, U., and Penzes, P. (2012) [Control of interneuron dendritic growth through NRG1/erbB4-mediated kalirin-7 disinhibition](#). *Mol Psychiatry* **17**, 1, 99-107
14. Cahill, M. E., Remmers, C., Jones, K. A., Xie, Z., Sweet, R. A., and Penzes, P. (2013) [Neuregulin1 signaling promotes dendritic spine growth through kalirin](#). *J Neurochem* **126**, 625-635
15. Jones, K. A., Srivastava, D. P., Allen, J. A., Strachan, R. T., Roth, B. L., and Penzes, P. (2009) [Rapid modulation of spine morphology by the 5-HT2A serotonin receptor through kalirin-7 signaling](#). *Proc Natl Acad Sci U S A* **106**, 19575-19580
16. Cahill, M. E., Xie, Z., Day, M., Photowala, H., Barbolina, M. V., Miller, C. A., Weiss, C., Radulovic, J., Sweatt, J. D., Disterhoft, J. F., Surmeier, D. J., and Penzes, P. (2009) [Kalirin regulates cortical spine morphogenesis and disease-related behavioral phenotypes](#). *Proc Natl Acad Sci U S A* **106**, 13058-13063
17. Ma, X. M., Kiraly, D. D., Gaier, E. D., Wang, Y., Kim, E. J., Levine, E. S., Eipper, B. A., and Mains, R. E. (2008) [Kalirin-7 is required for synaptic structure and function](#). *J Neurosci* **28**, 12368-12382
18. Youn, H., Jeoung, M., Koo, Y., Ji, H., Markesbery, W. R., Ji, I., and Ji, T. H. (2007) [Kalirin is under-expressed in Alzheimer's disease hippocampus](#). *J Alzheimers Dis* **11**, 385-397
19. Murray, P. S., Kirkwood, C. M., Gray, M. C., Ikonomic, M. D., Paljug, W. R., Abrahamson, E. E., Henteleff, R. A., Hamilton, R. L., Kofler, J. K., Klunk, W. E., Lopez, O. L., Penzes, P., and Sweet, R. A. (2012) [beta-Amyloid 42/40 ratio and kalirin expression in Alzheimer disease with psychosis](#). *Neurobiol Aging* **33**, 2807-2816
20. Ratovitski, E. A., Alam, M. R., Quick, R. A., McMillan, A., Bao, C., Kozlovsky, C., Hand, T. A., Johnson, R. C., Mains, R. E., Eipper, B. A., and Lowenstein, C. J. (1999) [Kalirin inhibition of inducible nitric-oxide synthase](#). *J Biol Chem* **274**, 993-999
21. Zhao, L., Ma, Q. L., Calon, F., Harris-White, M. E., Yang, F., Lim, G. P., Morihara, T., Ubeda, O. J., Ambegaokar, S., Hansen, J. E., Weisbart, R. H., Teter, B., Frautschy, S. A., and Cole, G. M. (2006) [Role](#)

- of p21-activated kinase pathway defects in the cognitive deficits of Alzheimer disease. *Nat Neurosci* **9**, 234-242
22. Penzes, P., Beeser, A., Chernoff, J., Schiller, M. R., Eipper, B. A., Mains, R. E., and Huganir, R. L. (2003) Rapid induction of dendritic spine morphogenesis by trans-synaptic ephrinB-EphB receptor activation of the Rho-GEF kalirin. *Neuron* **37**, 263-274
  23. Lesch, K. P., Timmesfeld, N., Renner, T. J., Halperin, R., Roser, C., Nguyen, T. T., Craig, D. W., Romanos, J., Heine, M., Meyer, J., Freitag, C., Warnke, A., Romanos, M., Schafer, H., Walitza, S., Reif, A., Stephan, D. A., and Jacob, C. (2008) Molecular genetics of adult ADHD: converging evidence from genome-wide association and extended pedigree linkage studies. *J Neural Transm (Vienna)* **115**, 1573-1585
  24. Xie, Z., Photowala, H., Cahill, M. E., Srivastava, D. P., Woolfrey, K. M., Shum, C. Y., Huganir, R. L., and Penzes, P. (2008) Coordination of synaptic adhesion with dendritic spine remodeling by AF-6 and kalirin-7. *J Neurosci* **28**, 6079-6091
  25. Mains, R. E., Kiraly, D. D., Eipper-Mains, J. E., Ma, X. M., and Eipper, B. A. (2011) Kalrn promoter usage and isoform expression respond to chronic cocaine exposure. *BMC Neurosci* **12**, 20
  26. Kiraly, D. D., Ma, X. M., Mazzone, C. M., Xin, X., Mains, R. E., and Eipper, B. A. (2010) Behavioral and morphological responses to cocaine require kalirin7. *Biol Psychiatry* **68**, 249-255
  27. Colomer, V., Engelender, S., Sharp, A. H., Duan, K., Cooper, J. K., Lanahan, A., Lyford, G., Worley, P., and Ross, C. A. (1997) Huntingtin-associated protein 1 (HAP1) binds to a Trio-like polypeptide, with a rac1 guanine nucleotide exchange factor domain. *Hum Mol Genet* **6**, 1519-1525
  28. Tsai, Y. C., Riess, O., Soehn, A. S., and Nguyen, H. P. (2012) The Guanine nucleotide exchange factor kalirin-7 is a novel synphilin-1 interacting protein and modifies synphilin-1 aggregate transport and formation. *PLoS One* **7**, e51999
  29. Li, W., Li, Q. J., and An, S. C. (2010) Preventive effect of estrogen on depression-like behavior induced by chronic restraint stress. *Neurosci Bull* **26**, 140-146
  30. Krug, T., Manso, H., Gouveia, L., Sobral, J., Xavier, J. M., Albergaria, I., Gaspar, G., Correia, M., Viana-Baptista, M., Simoes, R. M., Pinto, A. N., Taipa, R., Ferreira, C., Fontes, J. R., Silva, M. R., Gabriel, J. P., Matos, I., Lopes, G., Ferro, J. M., Vicente, A. M., and Oliveira, S. A. (2010) Kalirin: a novel genetic risk factor for ischemic stroke. *Hum Genet* **127**, 513-523
  31. Beresewicz, M., Kowalczyk, J. E., and Zablocka, B. (2008) Kalirin-7, a protein enriched in postsynaptic density, is involved in ischemic signal transduction. *Neurochem Res* **33**, 1789-1794

**We respectfully request that this document is cited using the DOI value as given above if the content is used in your work.**